

# HAS3 Rabbit pAb

货号: B13128

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  | <b>WB:</b> Homo sapiens                                                                                |
| 应用    | <a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF/ICC</a>                                          |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200                    |
| 理论分子量 | 31kDa/62kDa                                                                                            |
| 实测分子量 | 63KDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | PC-3,Rat lung                                                                                          |
| 细胞定位  | Membrane,Multi-pass membrane protein                                                                   |
| 纯化    | Affinity purification                                                                                  |


抗原信息

|      |                                                                                                                                                                                                                                 |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 67-281 of human HAS3 (NP_619515.1).                                                                                                               |
| 序列   | HRRMRAGQALKLPSPRRGSVLALCIAAYQEDPDYLRKCLRSAQRISFPDLKVVMMVVDGNRQEDAYMLDIFHEVLGGTEQA<br>GFFVWRSNFHEAGEGETEASLQEGMDRVRDVVRASTFSCIMQKWGGKREVMYAFKALGDSVDYIQVCDSDTVLDPA<br>CTIEMLRVLEEDPQVGGVGGDVQPPGKGMAVEDDQVQAAQVRATEAWSVHQRHVSREQ |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyaluronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene family. Compared to the proteins encoded by other members of this gene family, this protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multiple transcript variants. |
| 基因ID  | 3038                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 基因名   | HAS3                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Swiss | O00219                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| 别名    | HAS3                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)