

DET1 Rabbit pAb

货号: B13170

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	63kDa/65kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,Mouse thymus
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 15-240 of human DET1 (NP_060466.2).
序列	HVSTIKPRRIQNQNVIHRLERRRISGKAGTHWHQVRVFHQNVFPNFTVVNVEKPPCFLRKFS PDGRYFI AFSSDQTSLEIY EYQGCQAAEDLLQGYEGEILSNGNDQRSVNIRGRLFERFFVLLHITNVAANGEHLNRECSLFTDDCRCVIVGSAAYLPDEP HPPFFEYVRNSESVTPNPRSPLEDYSLHIIDLHTGRLCDTRTFKCDKVVL SHNQGLYLYKNIL

靶点信息

研究背景	Component of the E3 ubiquitin ligase DCX DET1-COP1 complex, which is required for ubiquitination and subsequent degradation of target proteins. The complex is involved in JUN ubiquitination and degradation.
基因ID	55070

基因名	DET1
Swiss	Q7L5Y6
别名	DET1

产品验证



Western blot analysis of DET1 expressed in BT-474,Mouse thymus using DET1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn