

G0S2 Rabbit pAb

货号: B13171

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	11kDa
实测分子量	11kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse kidney
细胞定位	Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human G0S2 (NP_056529.1).
序列	METVQELIPLAKEMMAQKRKGKMKVLYVLGSVLALFGVVLGLMETVCSPFTAARRLRDQEAAVAELQAALERQALQKQAL QEKGKQQDTVLGGRALSNRQ

靶点信息

研究背景	Promotes apoptosis by binding to BCL2, hence preventing the formation of protective BCL2-BAX heterodimers.
基因ID	50486

基因名	G0S2
Swiss	P27469
别名	G0S2

产品验证



Western blot analysis of G0S2 expressed in Mouse lung,Mouse kidney using G0S2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn