

RAB11FIP1 Rabbit pAb

货号: B13196

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IP: Homo sapiens
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	15kDa/47kDa/65kDa/70kDa/137kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,COS-7
细胞定位	Cytoplasmic vesicle,Recycling endosome,phagosome membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 330-540 of human RAB11FIP1 (NP_079427.4).
序列	EAKGEIKDSSPSSSPKGFRRKKHLFSSTENLAAGSWKEPAEGGGLSSDRQLSESSTKDSLKSMTLPSYRPAPLVSGDLRENMAPANSEATKEAKESKKPESRRSSLLSLMTGKKDVAKGSEGENPLTVPGREKGMGMGVKPGEDASGPAEDLVRSEKDTAAVVSRRQGSSSLNLFEDVQITEPEAEPESEPRPPISSPRAPQTRAVKP

靶点信息

研究背景	This gene encodes one of the Rab11-family interacting proteins (Rab11-FIPs), which play a role in the Rab-11 mediated recycling of vesicles. The encoded protein may be involved in endocytic sorting, trafficking of proteins including integrin subunits and epidermal growth factor receptor (EGFR), and transport between the recycling endosome and the trans-Golgi network. Alternative splicing results in multiple transcript variants. A pseudogene is described on the X chromosome.
------	--

基因ID	80223
基因名	RAB11FIP1
Swiss	Q6WKZ4
别名	RAB11FIP1;NOEL1A;RCP;rab11-FIP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn