

MYO18A Rabbit pAb

货号: B13204

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Gallus gallus
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	180kDa/196kDa/226kDa/231kDa/233kDa
实测分子量	260kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,Jurkat,Mouse brain,Mouse liver,Rat brain
细胞定位	Cytoplasm,Endoplasmic reticulum-Golgi intermediate compartment,Golgi apparatus,cytoskeleton,trans-Golgi network
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1970-2054 of human MYO18A (NP_510880.2).
序列	SDVDSELED RVDGVKSWLSKNKGPSKAASDDGSLKSSSPTS YWKS LAPDRSDDEHDPLDNTSRPRYSHSYLS DSDTEAKLTETNA

靶点信息

研究背景	The protein encoded by this gene can bind GOLPH3, linking the Golgi to the cytoskeleton and influencing Golgi membrane trafficking. The encoded protein is also part of a complex that assembles lamellar actomyosin bundles and may be required for cell migration.
基因ID	399687

基因名	MYO18A
Swiss	Q92614
别名	MYO18A;MAJN;MYSPDZ;SP-R210;SPR210

产品验证



Western blot analysis of MYO18A expressed in Mouse brain,Rat brain,U-251MG using MYO18A Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn