

BRF2 Rabbit pAb

货号: B13206

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IHC: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:100
理论分子量	44kDa/46kDa
实测分子量	62KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C6,Mouse testis,Rat lung
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 180-419 of human BRF2 (NP_060780.2).
序列	FKLFQASPSVPAKYVEDKEKMLSRTMQLVELANETWLV TGRHLPVITAATFLAWQSLQPADRLSCSLARFCKLANVDLPY PASSRLQELLAVLLRMAEQLAWLRVLR LDKRSVVKHIGDLLQHRQSLVRSAFRDGTAEVETREKEPPGWGGQGEGEVG NNSLGLPQGKRPA SPALLLPPCMLKSPKRICVPPVSTVTGDENISDSEIEQYL RTPQEV RDFQRAQAARQAATSVPNPP

靶点信息

研究背景	This gene encodes one of the multiple subunits of the RNA polymerase III transcription factor complex required for transcription of genes with promoter elements upstream of the initiation site. The product of this gene, a TFIIIB-like factor, is directly recruited to the TATA-box of polymerase III small nuclear RNA gene promoters through its interaction with the TATA-binding protein.
------	---

基因ID	55290
基因名	BRF2
Swiss	Q9HAW0
别名	BRF2;BRFU;TFIIIB50

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn