

MTERFD3 Rabbit pAb

货号: B13220

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	44kDa
实测分子量	44KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse heart,Rat heart
细胞定位	Mitochondrion,Mitochondrion matrix,mitochondrion nucleoid
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human MTERF D3 (NP_079474.2).
序列	MLWKLLLRSQSCRLCSFRKMRSPPKYRPFLACFTYTTDKQSSKENTRTVEKLYKCSVDIRKIRRLKGWVLLLEDETYVEEIANI LQELGADETAVASILERCPEAIVCSPTAVNTQRKLWQLVCKNEEELIKLIEQFPESFFTIKDQENQKLNQVQFFQELGLKNVVIS RLLTAAPNVFHNPPVEKNKQMVRLQESYLDVGGSEANMKVWLLKLLSQNPFIILLNSPTAIKETLEFLQEQQGFTSFEILQLLSK

靶点信息

研究背景	Binds mitochondrial DNA and plays a role in the regulation of transcription of mitochondrial mRNA and rR NA species.
基因ID	80298

基因名	MTERF2
Swiss	Q49AM1
别名	MTERF2;MTERFD3;mTERFL

产品验证



Western blot analysis of MTERFD3 expressed in Mouse brain,Mouse heart,Rat heart using MTERFD3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn