

TADA2A Rabbit pAb

货号: B13229

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	36kDa/51kDa
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,Mouse brain,Mouse skeletal muscle,Rat brain
细胞定位	Chromosome,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-443 of human TADA2 A (NP_001159577.2).
序列	MDRLGSFSNDPSDKPPCRGCSSYLMEPYIKCAECGPPPPFLCLQCFTRGFEYKKHQSDHTYEIMTSDFPVLDPSTWTAQEE MALLEAVMDCGFGNWQDVANQMCTKTKEECEKHYMKHFINNPLFASTLLNLKQAEAAKTADTAIPFHSTDDPPRPTFDSL LSRDMAGYMPARADFIEEFDNYAEWDLRDIDFVEDDSDILHALKMAVVDIYHSRLKERQRRKKIIRDHGLINLRKFQLMER RYPKEVQDLYETMRRFARIVGPVEHDKFIESHALEFELRREIKRLQEYRTAGITNFCSTARTYDHLKKTREEERLKRTMLSEVLQ YIQDSSACQQLRRQADIDSGLSIPMASNSGRRSAPPLNLTGLPGTEKLEKEKELCQMVRLVPGAYLEYKSALLNECN KQGGRLAQARALIKIDVKNKTRKIYDFLIREGYITKG

靶点信息

研究背景	Many DNA-binding transcriptional activator proteins enhance the initiation rate of RNA polymerase II-mediated gene transcription by interacting functionally with the general transcription machinery bound at the basal promoter. Adaptor proteins are usually required for this activation, possibly to acetylate and destabilize nucleosomes, thereby relieving chromatin constraints at the promoter. The protein encoded by this gene is a transcriptional activator adaptor and has been found to be part of the PCAF histone acetylase complex. Several alternatively spliced transcript variants encoding different isoforms of this gene have been described, but the full-length nature of some of these variants has not been determined.
基因ID	6871
基因名	TADA2A
Swiss	O75478
别名	TADA2A;ADA2;ADA2A;KL04P;TADA2L;hADA2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn