

# SLCO1A2 Rabbit pAb

货号: B13230

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>WB:</b> Rattus norvegicus                                                                         |
| 应用    | <a href="#">WB</a>                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:1000 - 1:5000                                                                           |
| 理论分子量 | 64kDa/74kDa                                                                                          |
| 实测分子量 | 74KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Rat liver                                                                                            |
| 细胞定位  | Cell membrane,Multi-pass membrane protein                                                            |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                        |
|------|------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 606-670 of human SLCO 1A2 (NP_602307.1). |
| 序列   | AALRGSSFVPALIIILLRKCHLPGENASSGTELIETKVKGKENECKDIYQKSTVLKDDDELKTKL                                                      |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                    |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a sodium-independent transporter which mediates cellular uptake of organic ions in the liver. Its substrates include bile acids, bromosulphophthalein, and some steroidal compounds. The protein is a member of the SLC21A family of solute carriers. Alternative splicing of this gene results in multiple transcript variants. |
| 基因ID | 6579                                                                                                                                                                                                                                                                                                                                               |

|       |                                     |
|-------|-------------------------------------|
| 基因名   | SLC01A2                             |
| Swiss | P46721                              |
| 别名    | SLC01A2;OATP;OATP-A;OATP1A2;SLC21A3 |

产品验证



Western blot analysis of SLC01A2 expressed in Rat liver using SLC01A2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)