

MARK1 Rabbit pAb

货号: B13232

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	72kDa/84kDa/89kDa
实测分子量	95kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U-87MG
细胞定位	Cell membrane,Cytoplasm,Peripheral membrane protein,cytoskeleton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 560-700 of human MARK1 (NP_061120.3).
序列	HPIKVTLP TI KDGSEAYRPGTTQ RV PAASPSAHSISTATPD RT FRPGSSSRSTFHGEQLRERRSVAYNGPPASPSHETGAF AHARRGTSTGIISKITSKFVRRDPSEGEASGR TD TSRSTSGEPKERDKEEGKDSK PR SL

靶点信息

研究背景	Serine/threonine-protein kinase. Involved in cell polarity and microtubule dynamics regulation. Phosphorylates DCX, MAP2 and MAP4. Phosphorylates the microtubule-associated protein MAPT/TAU. Involved in cell polarity by phosphorylating the microtubule-associated proteins MAP2, MAP4 and MAPT/TAU at KXGS motifs, causing detachment from microtubules, and their disassembly. Involved in the regulation of neuronal migration through its dual activities in regulating cellular polarity and microtubule dynamics, possibly by phosphorylating and regulating DCX. Also acts as a positive regulator of the Wnt signaling pathway, probably by mediating phosphorylation of dishevelled proteins (DVL1, DVL2 and/or DVL3).
基因ID	4139
基因名	MARK1
Swiss	Q9P0L2
别名	MARK1;MARK;Par-1c;Par1c

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn