

RNF11 Rabbit pAb

货号: B13249

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus ELISA: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	17kDa
实测分子量	27kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Rat liver
细胞定位	Cytoplasm,Early endosome,Nucleus,Recycling endosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-154 of human RNF11 (NP_055187.1).
序列	MGNCLKSPTSDDISLLHESQSDRASFGEGTEPDQEPPPPYQEQVPVPVYHPTPSQTRLATQLTEEEQIRIAQRIGLIQHLPK GVYDPGRDGSEKKIRECVICMMDFVYGDPIRFLPCMHIYHLCIDDWLMRSFTCPSCMEPVDAALLSSYETN

靶点信息

研究背景	The protein encoded by this gene contains a RING-H2 finger motif, which is known to be important for protein-protein interactions. The expression of this gene has been shown to be induced by mutant RET proteins (MEN2A/MEN2B). The germline mutations in RET gene are known to be responsible for the development of multiple endocrine neoplasia (MEN).
基因ID	26994

基因名	RNF11
Swiss	Q9Y3C5
别名	RNF11;CGI-123;SID1669

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn