

USF2 Rabbit pAb

货号: B13259

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	ChIP: Homo sapiens WB: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	23kDa/30kDa/36kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 101-200 of human USF2 (NP_003358.1).
序列	VVSTAAFAGGQQAVTQVGVDGAAQRPGPAAASVPPGPAAPFLAVIQNPFSNGGSPAAEAVSGEARFAYFPASSVGDTT AVSVQTTDQSLQAGGQFYVMM

靶点信息

研究背景	This gene encodes a member of the basic helix-loop-helix leucine zipper family of transcription factors. The encoded protein can activate transcription through pyrimidine-rich initiator (Inr) elements and E-box motifs and is involved in regulating multiple cellular processes.
基因ID	7392

基因名	USF2
Swiss	Q15853
别名	USF2;FIP;bHLHb12

产品验证



Western blot analysis of USF2 expressed in Rat liver using USF2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn