

CEACAM7 Rabbit pAb

货号: B13263

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens RT-qPCR: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa/29kDa
实测分子量	29kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Rat liver
细胞定位	Cell membrane,GPI-anchor,Lipid-anchor
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 36-242 of human CEACAM7 (NP_008821.2).
序列	TNIDVVPFNVAEGKEVLLVVHNESQNLGYGNWYKGERVHANYRIIGYVKNISQENAPGPAHNGRETIYPNGTLLIQNVTHND AGIYTLHVIKENLVNEEVTRQFYVFSEPPKPSITSNNFNPVENKDIVLTCQPETQNTTYLWWVNNQSLLVSPRLLLSTDNRT LVLLSATKNDIGPYECEIQNPVGASRSDPVTNLNVRYESVQAS

靶点信息

研究背景	This gene encodes a cell surface glycoprotein and member of the carcinoembryonic antigen (CEA) family of proteins. Expression of this gene may be downregulated in colon and rectal cancer. Additionally, lower expression levels of this gene may be predictive of rectal cancer recurrence. This gene is present in a CEA family gene cluster on chromosome 19. Alternative splicing results in multiple transcript variants.
------	---

基因ID	1087
基因名	CEACAM7
Swiss	Q14002
别名	CEACAM7;CGM2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn