

# RAB5B Rabbit pAb

货号: **B13311**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	<a href="#">WB</a> <a href="#">IF/ICC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	19kDa/23kDa
实测分子量	24kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SKOV3,Mouse ovary,Rat brain
细胞定位	Cell membrane,Cytoplasmic side,Early endosome membrane,Lipid-anchor,Melanosome
纯化	Affinity purification

## 抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RAB5B (NP_002859.1).
序列	YDITNQETFARAKTWVKELQRQASPSIVIALAGNKADLANKRMVEYEEAQAYADDNSLLFMETSAKTAMNVNDLFLAIKKLPKSEPQNLGGAAGRSRGVD

## 靶点信息

研究背景	Protein transport. Probably involved in vesicular traffic.
基因ID	5869
基因名	RAB5B

Swiss	P61020
别名	RAB5B

## 产品验证



Western blot analysis of RAB5B expressed in Mouse ovary,Rat brain,HeLa using RAB5B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)