

# SETDB2 Rabbit pAb

货号: B13321

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Mouse,Rat                                                                                            |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>WB:</b> Homo sapiens                                                                              |
| 应用    | <a href="#">WB</a>                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000                                                                            |
| 理论分子量 | 80kDa/81kDa                                                                                          |
| 实测分子量 | 100kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse heart,Rat heart                                                                                |
| 细胞定位  | Chromosome,Nucleus                                                                                   |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                        |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 148-378 of human SETD B2 (NP_114121.2).                                                                                                                  |
| 序列   | SGACLMKMPLNLKGENPLQLPIKCHFQRRHAKTNSHSSALHVSykTPCGRSLRNVEEVFRYLLETECNFLFTDNFSFNTYVQLARNYPKQKEVVSDVDISNGVESVPISFCNEIDSRKLPQFKYRKTVWPRAYNLTNFSSMFTDSCDCSEGCIDITKCACLQLTARNAKTSPLSSDKITTYKYKRLQRQIPTGIYECSSLCKCNRQLCQNRVVQHGPQVRLQVFKTEQKG |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of a family of proteins that contain a methyl-CpG-binding domain (MBD) and a SET domain and function as histone methyltransferases. This protein is recruited to heterochromatin and plays a role in the regulation of chromosome segregation. This region is commonly deleted in chronic lymphocytic leukemia. Naturally-occurring readthrough transcription occurs from this gene to the downstream PHF11 (PHD finger protein 11) gene. Alternative splicing results in multiple transcript variants. |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                  |
|-------|----------------------------------|
| 基因ID  | 83852                            |
| 基因名   | SETDB2                           |
| Swiss | Q96T68                           |
| 别名    | SETDB2;C13orf4;CLLD8;CLLL8;KMT1F |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)