

NDP52 Rabbit pAb

货号: B13324

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Saccharomyces cerevisiae
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	43kDa/47kDa/52kDa/54kDa/55kDa
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse pancreas
细胞定位	Cytoplasm,Cytoplasmic vesicle,Peripheral membrane protein,autophagosome membrane,cytoskeleton,pe rinuclear region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 137-446 of human NDP 52 (NP_005822.1).
序列	TTQGEVEEIEQHNLCKENQELKDSCISLQKQNSDMQAELQKKQELETLSINKKLELKVKEQKDYWETELLQLKEQNQ KMSSENEKMGIRVDQLQAQLSTQEKEMEKLVGDDQDKTEQLEQLKKENDHFLSLTEQRKDQKKLEQTVEQMKQNETT AMKKQQELMDENFDLSKRLSENEIICNALQRQKERLEGENDLLKRENSRLLSYMGLDFNSLPYQVPTSDEGGARQNPGLA YGNPYSGIQESSSPSLSIKKCPICKADDICDHTLEQQMQPLCFNCPICDKIFPATEKQIFEDHVFCHSL

靶点信息

研究背景	This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as a receptor for ubiquitin-coated bacteria and plays an important role in innate immunity by mediating macroautophagy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	---

基因ID	10241
基因名	CALCOCO2
Swiss	Q13137
别名	CALCOCO2;NDP52

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn