

RNF8 Rabbit pAb

货号: B13333

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	10kDa/50kDa/55kDa
实测分子量	60KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,Mouse brain,Mouse liver,Rat brain
细胞定位	Chromosome,Midbody,Nucleus,telomere
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RNF8 (NP_003949.1).
序列	MGEPGFFVTGDRAGGRSWCLRRVGM SAGWLLLEDGCEVTVGRGFGVTYQLVSKICPLMISRNHCVLKQNPEGQWTIM DNKSLNGVWLNRRARLEPLRVYSIHQGDYIQLGVPLENKENAEYEVTEEDWETIYPCLSPKNDQMIEKNKELRTKRKFSLD ELAGPGAEGPSNLKSKINKVSCESGQPVKSQGKGEVASTPSDNLDPKLTALEPSKTTGAPIYPGFPKVT EVHHHEQKASNSS ASQRSLQMFKVTMSRILRLK

靶点信息

研究背景	The protein encoded by this gene contains a RING finger motif and an FHA domain. This protein has been shown to interact with several class II ubiquitin-conjugating enzymes (E2), including UBE2E1/UBCH6, UBE 2E2, and UBE2E3, and may act as an ubiquitin ligase (E3) in the ubiquitination of certain nuclear proteins. This protein is also known to play a role in the DNA damage response and depletion of this protein causes cell growth inhibition and cell cycle arrest. Alternative splicing results in multiple transcript variants.
------	--

基因ID	9025
基因名	RNF8
Swiss	O76064
别名	RNF8;hRNF8

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn