

ING5 Rabbit pAb

货号: B13337

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:100
理论分子量	25kDa/27kDa
实测分子量	30KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ING5 (NP_115705.2).
序列	MATAMYLEHYLDSIENLPCELQRNFQLMRELDQRTEDKKAIEDILAAEYISTVKTLSPDQRVERLQKIQNAYSKCKEYSDDK VQLAMQTYEMVDKHIRRLDADLARFEADLKDKMEGSDFESSGGRLKKGRGQKEKRGSRGRGRRTSEEDTPKKKKHK GGSEFTDILSVHPSDVLDMPVDPNEPTYCLCHQVSYGEMIGCDNPDCPIEFWFHFACVDLT

靶点信息

研究背景	This gene encodes a tumor suppressor protein that inhibits cell growth and induces apoptosis. This protein contains a PHD-type zinc finger. It interacts with tumor suppressor p53 and p300, a component of the histone acetyl transferase complex, suggesting a role in transcriptional regulation. Alternative splicing and the use of multiple promoters and 3' ends results in multiple transcript variants.
------	--

基因ID	84289
基因名	ING5
Swiss	Q8WYH8
别名	ING5;p28ING5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn