

Sall4 Rabbit pAb

货号: B13352

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	65kDa/112kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,Mouse heart
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human Sall4 (NP_065169.1).
序列	MSRRKQAKPQHINSEEDQGEQQPQQQTPEFADAAPAAPAAGELGAPVNHHPGNDEVASEDEATVKRLRREETHVCEKCCAEFFSISEFLEHKKNCTKNPPVLIMNDSEGPVPSDFSGAVLSHQPTSPGSKDCHRENGGSSEDMKEKPD AESVVYLKTE TALPPTPQDISYLAKGKVANTNVTLQALRGTKVAVNQRSADALPAPVPGANSIPWVLEQILC

靶点信息

研究背景	This gene encodes a zinc finger transcription factor thought to play a role in the development of abducens motor neurons. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS). Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	57167

基因名	SALL4
Swiss	Q9UJQ4
别名	SALL4;DRRS;HSAL4;ZNF797

产品验证



Western blot analysis of Sall4 expressed in Jurkat,HeLa,Mouse heart using Sall4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn