

FKBP38/FKBP8 Rabbit pAb

货号: B13359

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100 IP: 1:50 - 1:200
理论分子量	27kDa/44kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,THP-1,BT-474,HeLa,HL-60,Mouse liver,Mouse brain,Mouse eye,Rat brain
细胞定位	Cytoplasmic side,Mitochondrion membrane,Single-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-270 of human FKBP38/FKBP38/FKBP8 (NP_036313.3).
序列	QPPAEAEQPGALAREFLAAMEPEPAPAPAPEEWLDILGNGLLRKKTLPVPGPPGSSRPVKGQVVTVHLQTSLENGTRVQEE PELVFTLGDCDVIQALDLSVPLMDVGETAMVTADSKYCYGPQGSRSPIPPHAALCLEVTLKTAVDGPDLEMLTGQERVAL ANRKRECGNAHYQRADFVLAANSYDLAIKAITSSAKVDMTFEEEAQLLQ

靶点信息

研究背景	The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. Unlike the other members of the family, this encoded protein does not seem to have PPlase/rotamase activity. It may have a role in neurons associated with memory function.
------	--

基因ID	23770
基因名	FKBP8
Swiss	Q14318
别名	FKBP8;FKBP38;FKBPr38

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn