

TOX Rabbit pAb

货号: B13363

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Mus musculus
应用	WB IF/ICC IHC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IHC: 1:50 - 1:200
理论分子量	57kDa
实测分子量	57kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 297-526 of human TOX (NP_055544.1).
序列	SMWDGLGEEQKQVYKKKTEAAKKEYLKQLAAYRASLVSKSYSEPVDVKTSQPPQLINSKPSVFHGPSQAHSALYLSSHYY QQPGMNPHLTAMHPSLPRNIAPKPNQMPVTVSIANMAVSPPPPLQISPLHQHLNMQQHQPLTMQQPLGNQLPMQVQ SALHSPTMQQGFTLQPDYQTIINPTSTAAQVVTQAMEYVRSGRNPPPQPVVDWNNDYCSSGGMQRDKALYLT

靶点信息

研究背景	The protein encoded by this gene contains a HMG box DNA binding domain. HMG boxes are found in many eukaryotic proteins involved in chromatin assembly, transcription and replication. This protein may function to regulate T-cell development.
------	--

基因ID	9760
基因名	TOX
Swiss	O94900
别名	TOX;TOX1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn