

MUS81 Rabbit pAb

货号: B13387

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	61kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat
细胞定位	Nucleus,nucleolus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human MUS81 (NP_079404.3).
序列	MAAPVRLGRKRPLPACPNPLFVRWLTEWRDEATRSRRRTRFVFQKALRSLRRYPLPLRSGKEAKILQHFGDGLCRMLDERLQRHRTSGGDHAPDSPGENSPAPQGRLAEVQDSSMPVPAQPKAGGSGSYWPARHSGARVILLVLYREHLNPNGHHFLTKEELLQRCAQKSPRVAPGSARPWPALRSLLHRNLVLRTHQPARYSLTPEGLELAQKLAESEGLSLLNVGIGPKEPPGEETAVPGAASAELASEAGVQQQPLELRPGEYRVLLCVDIGETRGGGHRPELLRELQRLHVTHTVRKLHVGDVWVAQETNPRDPA NPGELVLDHIVERKRLDDLCSSIIDGRFR

靶点信息

研究背景	This gene encodes a structure-specific endonuclease which belongs to the XPF/MUS81 endonuclease family and plays a critical role in the resolution of recombination intermediates during DNA repair after inter-strand cross-links, replication fork collapse, and DNA double-strand breaks. The encoded protein associates with one of two closely related essential meiotic endonuclease proteins (EME1 or EME2) to form a complex that processes DNA secondary structures. It contains an N-terminal DEAH helicase domain, an excision repair cross complementation group 4 (ERCC4) endonuclease domain, and two tandem C-terminal helix-hairpin-helix domains. Mice with a homozygous knockout of the orthologous gene have significant meiotic defects including the failure to repair a subset of DNA double strand breaks.
基因ID	80198
基因名	MUS81
Swiss	Q96NY9
别名	MUS81;SLX3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn