

TPPP3 Rabbit pAb

货号: B13390

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC, WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	18kDa
实测分子量	19KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,A-549,Mouse spinal cord,Rat lung
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human TPPP3 (NP_057224.2).
序列	MAASTDMAGLEESFRKFAIHGDPKASGQEMNGKNWAKLCKDCKVADGKSVTGTDVDIVFSKVKGKSARVINYEFEKKAL EELATKRFGKGSKEEAFDAICQLVAGKEPANVGVTAKTGGAVDRLTDSRYTGSHKERFDESGKGKGIAGRQDILDDSG YVSAYKNAGTYDAKVKK

靶点信息

研究背景	Regulator of microtubule dynamic that has microtubule bundling activity. Required for embryo implantati on; possibly by regulating beta-catenin (By similarity. Also required for decidualization via regulation of b eta-catenin.
基因ID	51673

基因名	TPPP3
Swiss	Q9BW30
别名	TPPP3;CGI-38;TPPP/p20;p20;p25gamma

产品验证



Western blot analysis of TPPP3 expressed in Mouse spinal cord,Rat lung,HT-29 using TPPP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn