

SLC26A3 Rabbit pAb

货号: B13393

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus,Sus scrofa
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	84kDa
实测分子量	85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H125
细胞定位	Apical cell membrane,Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 505-764 of human SLC26A3 (NP_000102.1).
序列	FPKCSTLANIGRTNIYKNKKDYDDMYEPEGVKIFRCPSPIYFANIGFFRRKLIDAVGFSPRLRLRKRNKALRKIRKLQKQGLLQVTPKGFICTVDTIKDSDEELDNNQIEVLDPINTTDLPFHIDWNDDLPLNIEVPKISLHSLILDFSAVSFLDVSSVRGLKSILQEFIRIKVDVYIVGTDDDFIEKLNRYEFFDGEVKSSIFFLTIHDAVLHILMKKDYSTSKFNPSQEKDGKIDFTINTNGGLRNRVYEV PVETKF

靶点信息

研究背景	The protein encoded by this gene is a transmembrane glycoprotein that transports chloride ions across the cell membrane in exchange for bicarbonate ions. It is localized to the mucosa of the lower intestinal tract, particularly to the apical membrane of columnar epithelium and some goblet cells. The protein is essential for intestinal chloride absorption, and mutations in this gene have been associated with congenital chloride diarrhea.
------	--

基因ID	1811
基因名	SLC26A3
Swiss	P40879
别名	SLC26A3;CLD;DRA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn