

OSMR Rabbit pAb

货号: B13401

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	39kDa/110kDa
实测分子量	111KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-330 of human OSMR (NP_003990.1).
序列	LPLTPVSLKVSTNSTRQSLHLQWTVHNLPHYQELKMFVQIQISRIETSNVIWVGNYSTTVKWNQVLHWSWESELPLECATHFVRIKSLVDDAKFPEPNFWSNWSSWEEVSVQDSTGQDILFVFPKDKLVEEGTNVTICYVSRNIQNNVSCYLEGKQIHGEQLDPHVTAFLNLSVPFIRNKGFTNIYCEASQGNVSEGMKGIVLFVSKVLEEPKDFSCETEDFKTLHCTWDPGTDALGWSKQPSQSYTLFESFSGEKKLCTHKNWCNWQITQDSQETYNFTLIAENYLKRKRSVNILFNLTHR

靶点信息

研究背景	This gene encodes a member of the type I cytokine receptor family. The encoded protein heterodimerizes with interleukin 6 signal transducer to form the type II oncostatin M receptor and with interleukin 31 receptor A to form the interleukin 31 receptor, and thus transduces oncostatin M and interleukin 31 induced signaling events. Mutations in this gene have been associated with familial primary localized cutaneous amyloidosis. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	9180
基因名	OSMR
Swiss	Q99650
别名	OSMR;IL-31R-beta;IL-31RB;OSMRB;PLCA1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn