

NLRP12 Rabbit pAb

货号: B13404

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	31kDa/101kDa/107kDa/113kDa/114kDa/120kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,LO2
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 862-1062 of human NLR P12 (NP_001264055.1).
序列	LWLKICRLTAAACDELASTLSVNQSLRELDLSLNELGDLGVLLLCCEGLRHPTCKLQTLRLGICRLGSAACEGLSVVLQANHN LRELDLSFNDLGDWGLWLLAEGLQHPACRLQKLWLDSCGLTAKACENLYFTLGINQTLTDLYLTNNALGDTGVRLLCKRL SHPGCKLRVLWLFMDLNKMTHSRLAALRVTKPYLDIGC

靶点信息

研究背景	This gene encodes a member of the CATERPILLER family of cytoplasmic proteins. The encoded protein, w hich contains an N-terminal pyrin domain, a NACHT domain, a NACHT-associated domain, and a C-termin us leucine-rich repeat region, functions as an attenuating factor of inflammation by suppressing inflamma tory responses in activated monocytes. Mutations in this gene cause familial cold autoinflammatory syndr ome type 2. Alternative splicing results in multiple transcript variants.
------	---

基因ID	91662
基因名	NLRP12
Swiss	P59046
别名	NLRP12;CLR19.3;FCAS2;NALP12;PAN6;PYPAF7;RNO;RNO2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn