

INTS6 Rabbit pAb

货号: B13409

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	12kDa/98kDa/100kDa
实测分子量	117KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse testis
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human INTS6 (NP_036273.1).
序列	MPILLFLIDTSASMNQRSHLGTTYLDTAKGAVETFMKLRRDPASRGDRYMLVTFEPPYAIKAGWKENHATFMNELKNLQ AEGTLTLGQSLRTAFDLLNLRLVTGIDNYGQGRNPFFLEPAIIITITDGSKLTTTSGVQDELHLPLNSPLPGSELTKEPFRW DQRLFALVLRPLPGTMSVESEQLTGVPLDDSAITPMCEVTGGRSYSVCSPRMLNQCLESLVQKVQSGVVINFEKAGPDPSP VEDGQPDISRPFGSQPWHSCHKLIYVRPNPKTGVPI

靶点信息

研究背景	DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. The protein encoded by this gene is a DEAD box protein that is part of a complex that interacts with the C-terminus of RNA polymerase II and is involved in 3' end processing of snRNAs. In addition, this gene is a candidate tumor suppressor and is located in the critical region of loss of heterozygosity (LOH). Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	26512
基因名	INTS6
Swiss	Q9UL03
别名	INTS6;DBI-1;DDX26;DDX26A;DICE1;HDB;INT6;NotchI2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn