

ETFDH Rabbit pAb

货号: B13419

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	62kDa/68kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse heart,Mouse kidney,Rat liver,Rat heart
细胞定位	Mitochondrion inner membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 318-617 of human ETFD H (NP_004444.2).
序列	LVALGLVVGLDYQNPYLSPFREFQRWKHHPSIRPTLEGGKRIAYGARALNEGGFQSIPKLTFPGLLIGCSPGFMNVPKIKG THTAMKSGILAAESIFNQLTSENLSKTIGLHVTEYEDNLKNSWVWKELYSVRNIRPSCHGVLGVYGGMIYTGIFYWILRGM EPWTLKHKGSDFERLKPAKDCTPIEYPKPDGQISFDLLSSVALSGTNHEHDQPAHLTLRDDSI PVNRNLSIYDGPEQRFCP AGVYEFVPVEQGDGFRLLQINAQNCVHCKTCDIKDPSQININWVWVEGGGGPAYNGM

靶点信息

研究背景	This gene encodes a component of the electron-transfer system in mitochondria and is essential for electron transfer from a number of mitochondrial flavin-containing dehydrogenases to the main respiratory chain. Mutations in this gene are associated with glutaric acidemia. Alternatively spliced transcript variants that encode distinct isoforms have been observed.
------	---

基因ID	2110
基因名	ETFDH
Swiss	Q16134
别名	ETFDH;ETFQO;MADD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn