

FGF19 Rabbit pAb

货号: B13421

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	24kDa
实测分子量	22KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-216 of human FGF19 (NP_005108.1).
序列	LAFSDAGPHVHYGWGDPIRLRHLYTSGPHGLSSCFLRIRADGVVDCARGQSAHSLEIKAVALRTVAIKGVHSVRYLCMG ADGKMQGLLQYSEEDCAFEIEIRPDGYNVYRSEKHRLPVLSLSSAKQRQLYKNRGFLPLSHFLPMLPMVPEEPEDLRGHLES DMFSSPLETDSMDPFGGLVTGLEAVRSPSFEK

靶点信息

研究背景	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes including embryonic development cell growth, morphogenesis, tissue repair, tumor growth and invasion. This growth factor is a high affinity, heparin dependent ligand for FGFR4. Expression of this gene was detected only in fetal but not adult brain tissue. Synergistic interaction of the chick homolog and Wnt-8c has been shown to be required for initiation of inner ear development.
------	--

基因ID	9965
基因名	FGF19
Swiss	O95750
别名	FGF19

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn