

# BLMH Rabbit pAb

货号: B13425

## 产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  | <b>WB:</b> Homo sapiens                                                                                |
| 应用    | <a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF/ICC</a>                                          |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:20 - 1:50                    |
| 理论分子量 | 52kDa                                                                                                  |
| 实测分子量 | 50kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | NCI-H460,Rat thymus                                                                                    |
| 细胞定位  | Cytoplasm                                                                                              |
| 纯化    | Affinity purification                                                                                  |


## 抗原信息

|      |                                                                                                                                                                                                                                                                                                                 |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human BLMH ( NP_000377.1).                                                                                                                                                                                               |
| 序列   | MSSSGLNSEKVAALIQKLNSDPQFVLAQNVGTTDLLDICLKRA TVQRAQHVFQHAVPQEGKPITNQSSGRCWIFSCLN VMRLPFMKKLNIEEFESQSYLFFWDKVERCYFFLSAFVDTAQRKEPEDGRLVQFLLMNPANDGGQWDMLVNIVEKYGVI PKKCFPESYTTEATRRMNDILNHKMREFCIRLRNLVHSGATKGEISATQDVMMEEIFRVVCICLG NPPETFTWEYRDKDKNY QKIGPITPLEFYREHVKPLFNMEDKICLVNDPRPQHKYNKLYTVEYLSNMVGG RKTLY |

## 靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | Bleomycin hydrolase (BMH) is a cytoplasmic cysteine peptidase that is highly conserved through evolution; however, the only known activity of the enzyme is metabolic inactivation of the glycopeptide bleomycin (BLM), an essential component of combination chemotherapy regimens for cancer. The protein contains the signature active site residues of the cysteine protease papain superfamily. |
| 基因ID  | 642                                                                                                                                                                                                                                                                                                                                                                                                  |
| 基因名   | BLMH                                                                                                                                                                                                                                                                                                                                                                                                 |
| Swiss | Q13867                                                                                                                                                                                                                                                                                                                                                                                               |
| 别名    | BLMH;BH;BMH                                                                                                                                                                                                                                                                                                                                                                                          |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)