

TRIM27 Rabbit pAb

货号: B13440

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	41kDa/58kDa
实测分子量	59kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,mouse testis
细胞定位	Cytoplasm,Early endosome,Nucleus,PML body
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 214-513 of human TRIM 27 (NP_006501.1).
序列	SINGAITQFSCNISHLSSLIAQLEEKQQQPTRELLQDIGDTLSRAERIRIPEPWITPPDLQEKIHIFAQKCLFLTESLKQFTEKM QSDMEKIQELREAQLYSVDVTLDPTAYPSLILSDNLRQVRYSYLQQDLPDNPERFNLFPCLVGSPCFIAGRHYWEVEVGD KAKWTIGVCEDSVCRKGGVTSAPQNGFWAVSLWYGKEYWALTSPMTALPLRTPLQRVGIFLDYDAGEVSFYNVTERCHT FTFSHATFCGPVRPYFSLSYSGGKSAAPLIICPMSGIDGFSGHVGNHGHSMETSP

靶点信息

研究背景	This gene encodes a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to the nuclear matrix. It interacts with the enhancer of polycomb protein and represses gene transcription. It is also thought to be involved in the differentiation of male germ cells. Fusion of the N-terminus of this protein with the truncated C-terminus of the RET gene product has been shown to result in production of the ret transforming protein.
基因ID	5987
基因名	TRIM27
Swiss	P14373
别名	TRIM27;RFP;RNF76

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn