

SERPINB9 Rabbit pAb

货号: B13441

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:10 - 1:100
理论分子量	42kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,SKOV3,Mouse lung,Mouse spleen,Rat spleen,Rat thymus
细胞定位	Cytoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-376 of human SERPIN B9 (NP_004146.1).
序列	METLSNASGTFAIRLLKILCQDNPSHNVFCSPVSISSALAMVLLGAKGNTATQMAQALSLNTEEDIHRAFQSLLTEV NKAGT QYLLRTANRLFGEKTCQLSTFKESCLQFYHAELKELSFIRAAEESRKHINTWVSKKTEGKIEELLPGSSIDAETRLVLVNAIYF KGKWNPEPFDETYTREMPFKINQEEQRPVQMMYQEATFKLAHVGEVRAQLLELPYARKELSLLVLLPDDGVELSTVEKSLTF EKLTAWTKPDCMKSTEVEVLLPKFKLQEDYDMESVLRHLGIVDAFQQGKADLSAMSAERDLCLSKFVHKSFVEVNEEGTE AAAASSCFVVAECCMESGPRFCADHPFLFFIRHNRANSILFCGRFSSP

靶点信息

研究背景	This gene encodes a member of the serine protease inhibitor family which are also known as serpins. The encoded protein belongs to a subfamily of intracellular serpins. This protein inhibits the activity of the effector molecule granzyme B. Overexpression of this protein may prevent cytotoxic T-lymphocytes from eliminating certain tumor cells. A pseudogene of this gene is found on chromosome 6.
基因ID	5272
基因名	SERPINB9
Swiss	P50453
别名	SERPINB9;CAP-3;CAP3;PI-9;PI9;serpin B9

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn