

SMARCD1 Rabbit pAb

货号: B13453

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	53kDa/58kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat brain,Mouse thymus,Mouse heart,Mouse testis,Mouse brain,Jurkat,SW480,BT-474,SW620
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 216-515 of human SMA RCD1 (NP_003067.3).
序列	LEDSALSKYDATKQKRKFSSFFKSLVIELDKDLYGPDNHLVEWHRTATTQETDGFQVKRPGDVNVRCTVLLMLDYQPPQF KLDPRLARLLGIHTQTRPVIIQALWQYIKTHKLQDPHEREFVICDKYLQQIFESQRMKFSEIPQRLHALLMPPEPIIIINHVISVDP NDQKKTACYDIDVEVDDTLKTQMNSFLLSTASQQEIATLDNKIHETIETINQLKTQREFMLSFARDPQGFINDWLQSQCRD LKTMTDVVGNPEEERRAEFYFQPWAEAVCRYFYSKVQRRQELEQALGIRNT

靶点信息

研究背景	The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by alterin g the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent c hromatin remodeling complex SNF/SWI and has sequence similarity to the yeast Swp73 protein. Two tran script variants encoding different isoforms have been found for this gene.
------	---

基因ID	6602
基因名	SMARCD1
Swiss	Q96GM5
别名	SMARCD1;BAF60A;CRACD1;Rsc6p

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn