

G3BP2 Rabbit pAb

货号: B13473

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Mus musculus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	50kDa/54kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 140-280 of human G3BP 2 (NP_987100.1).
序列	DSEPELDEESEDEVEEEEQEERQPSPEPVQENANSGYYEAHPVTNGIEEPLIESSHEPEPEPESETKTEELKPQVEEKNLEELE EKSTTPPPAEPVSLPQEPPKPRVEAKPEVQSQPPRVREQRPRERPGFPFPRGPRPGRGD

靶点信息

研究背景	Scaffold protein that plays an essential role in cytoplasmic stress granule formation which acts as a platform for antiviral signaling.
基因ID	9908

基因名	G3BP2
Swiss	Q9UN86
别名	G3BP2

产品验证



Western blot analysis of G3BP2 expressed in HL-60 using G3BP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn