

IFNGR1 Rabbit pAb

货号: **B13502**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Rattus norvegicus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	21kDa/54kDa
实测分子量	80-90KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse liver,Rat liver
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 270-489 of human IFNG R1 (NP_000407.1).
序列	IKKINPLKEKSIILPKSLISVVRSATLETKPESKYVSLITSYQPFSLEKEVVC EEPLSPATVPGMHTEDNPGKVEHTEELSSITEVV TTEENIPDVVPGSHLTPIERESSPLSSNQSEPGSIALNSYHSRNCSESDHSRNGFDTDSSCLES HSSLSDSEFPNNKGEI KTEGQELITVIKAPTSFGYDKPHVLVDLLVDDSGKESLIGYRPTEDSKEFS

靶点信息

研究背景	This gene (IFNGR1) encodes the ligand-binding chain (alpha) of the gamma interferon receptor. Human in terferon-gamma receptor is a heterodimer of IFNGR1 and IFNGR2. A genetic variation in IFNGR1 is associ ated with susceptibility to Helicobacter pylori infection. In addition, defects in IFNGR1 are a cause of men delian susceptibility to mycobacterial disease, also known as familial disseminated atypical mycobacterial infection.
------	--

基因ID	3459
基因名	IFNGR1
Swiss	P15260
别名	IFNGR1;CD119;IFNGR;IMD27A;IMD27B

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn