

CDKN2B Rabbit pAb

货号: B13507

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	8kDa/14kDa
实测分子量	15KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,THP-1(negative),Mouse small intestine
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 59-138 of human CDKN2B (NP_004927.2).
序列	ARVAELLLLHGAEPNCADPATLTRPVHDAAREGFDLTVVLHRAGARLDVRDAWGRLPVDLAEERGHDRVAGYLRATG D

靶点信息

研究背景	This gene lies adjacent to the tumor suppressor gene CDKN2A in a region that is frequently mutated and deleted in a wide variety of tumors. This gene encodes a cyclin-dependent kinase inhibitor, which forms a complex with CDK4 or CDK6, and prevents the activation of the CDK kinases, thus the encoded protein functions as a cell growth regulator that controls cell cycle G1 progression. The expression of this gene was found to be dramatically induced by TGF beta, which suggested its role in the TGF beta induced growth inhibition. Two alternatively spliced transcript variants of this gene, which encode distinct proteins, have been reported.
------	--

基因ID	1030
基因名	CDKN2B
Swiss	P42772
别名	CDKN2B;CDK4I;INK4B;MTS2;P15;TP15;p15INK4b

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn