

Hemopexin (HPX) Rabbit pAb

货号: B13520

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,HepG2,BT-474,Mouse liver,Mouse pancreas
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 24-245 of human Hemopexin (Hemopexin (HPX)) (NP_000604.1).
序列	TPLPPTSAHGNVAEGETKPDVDVTERCSDGWSFDATLDDNGTMLFFKGEFVWKSHKWDRELISERWKNFPSPVDAAF RQGHNSVFLIKGDKVWVYPPEKKEKGYPKLLQDEFFGIPSPLDAAVECHRGECQAEGVLFFQGDREWFWDLATGTMKER SWPAVGNCSSALRWLGRIYCFQGNQFLRFDVPVRGEVPPRYPRDVRDYFMPCPGRGHGHRNGTGHG

靶点信息

研究背景	This gene encodes a plasma glycoprotein that binds heme with high affinity. The encoded protein is an acute phase protein that transports heme from the plasma to the liver and may be involved in protecting cells from oxidative stress.
基因ID	3263

基因名	HPX
Swiss	P02790
别名	HPX;HX;hemopexin

产品验证



Western blot analysis of Hemopexin (HPX) expressed in Mouse liver,HL-60,HepG2 using Hemopexin (HPX) Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysate s/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Ki t. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn