

CTSA Rabbit pAb

货号: B13533

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens,Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	52kDa/54kDa
实测分子量	54KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse brain,Rat kidney
细胞定位	Lysosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 210-450 of human CTSA (NP_001121167.1).
序列	SYEQNDNSLVYFAYYHGLLGNRLWSSLQTHCCSQNKCIFYDNKDLECVTNLQEVARIVGNSGLNIYNLYAPCAGGVPSHF RYEKDTVVVQDLGNIFTRLPLKRMWHQALLRSGDKVRMDPPCTNTTAASTYLNPNPYVRKALNIPEQLPQWDMCNFLVNL QYRRLYRSMNSQYLKLLSSQKYQILLYNGDVDMACNFMGDEWFVDSLNLQKMEVQRRPWLVKYGDSGEQIAGFVKEFSH IAFL

靶点信息

研究背景	This gene encodes a member of the peptidase S10 family of serine carboxypeptidases. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate two chains that comprise the heterodimeric active enzyme. This enzyme possesses deamidase, esterase and carboxypeptidase activities and acts as a scaffold in the lysosomal multienzyme complex. Mutations in this gene are associated with galactosialidosis.
------	--

基因ID	5476
基因名	CTSA
Swiss	P10619
别名	CTSA;GLB2;GSL;NGBE;PPCA;PPGB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn