

SNRPE Rabbit pAb

货号: B13535

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	10kDa
实测分子量	11kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,HepG2,U-937
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human SNRPE (NP_003085.1).
序列	MAYRGQGQKVKVMVQPINLIFRYLQNRSRIQVWLYEQVMRIEGCIIGFDEYMNVLDDAEEIHSKTKSRKQLGRIMLKG DNITLLQSVSN

靶点信息

研究背景	The protein encoded by this gene is a core component of U small nuclear ribonucleoproteins, which are key components of the pre-mRNA processing spliceosome. The encoded protein plays a role in the 3' end processing of histone transcripts. This protein is one of the targets in the autoimmune disease systemic lupus erythematosus, and mutations in this gene have been associated with hypotrichosis. Several pseudogenes of this gene have been identified.
------	--

基因ID	6635
基因名	SNRPE
Swiss	P62304
别名	SNRPE;HYPT11;SME;Sm-E;snRNP-E

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn