

GNB5 Rabbit pAb

货号: B13586

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	31kDa/38kDa/43kDa
实测分子量	39kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HL-60,293T,Mouse brain,Rat brain
细胞定位	cytoplasm,cytosol,nucleus,presynapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human GNB5 (NP_006569.1).
序列	MATEGLHENETLASLKSEAESLKGKLEEEERAKLHDVELHQVAERVEALGQFVMKTRRTLKGHGKVL CMDWCKDKRRIVS SSQDGKVIVWDSFTTNKEHAVTMPCTWVMACAYAPSGCAIACGGLDNKC SVYPLTFDKNENMAAKKKS VAMHTNYLSA CSFTNSDMQILTASGDGTCALWDVESGQLLQSFHGHGADVLC

靶点信息

研究背景	Heterotrimeric guanine nucleotide-binding proteins (G proteins), which integrate signals between receptors and effector proteins, are composed of an alpha, a beta, and a gamma subunit. These subunits are encoded by families of related genes. This gene encodes a beta subunit. Beta subunits are important regulators of alpha subunits, as well as of certain signal transduction receptors and effectors. Alternatively spliced transcript variants encoding different isoforms exist.
------	---

基因ID	10681
基因名	GNB5
Swiss	O14775
别名	GNB5;GB5;IDDCA;LADCI

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn