

VARS Rabbit pAb

货号: B13592

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	33kDa/140kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,HL-60,K-562,SW620,Mouse brain,Mouse heart,Mouse liver,Rat liver
细胞定位	cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human VARS (NP_006286.1).
序列	MSTLYVSPHPDAFPSLRALIAARYGEAGEGPGWGGAHPRICLQPPPTSRTPFPPPRLPALQGGPGGLWVWGATAVAQLLW PAGLGGPGGSRAAVLVQQWVSYADTELIPAACGATLPALGLRSSAQDPQAVLGALGRALSPLEEWLRLHTYLAGEAPTLA DLAAVTALLLPFRYVLDPPARRIWNNVTRWVFVTCVRQPEFRAVLGEVVLVSGARPLSHQPGPEAPALPKTAAQLKKEAKKR EKLEKFQKQKIQQQQPPPGEEKPKPEKREKRDPGVITY

靶点信息

研究背景	Aminoacyl-tRNA synthetases catalyze the aminoacylation of tRNA by their cognate amino acid. Because of their central role in linking amino acids with nucleotide triplets contained in tRNAs, aminoacyl-tRNA synthetases are thought to be among the first proteins that appeared in evolution. The protein encoded by this gene belongs to class-I aminoacyl-tRNA synthetase family and is located in the class III region of the major histocompatibility complex.
基因ID	7407
基因名	VARS
Swiss	P26640
别名	VARS;G7A;VARS1;VARS2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn