

RALB Rabbit pAb

货号: B13597

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	23kDa/25kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,HepG2,Mouse lung,Rat lung
细胞定位	Cell membrane,Cytoplasmic side,Lipid-anchor
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human RALB (NP_002872.1).
序列	MAANKSKGQSSLALHKVIMVGSGGVGKSALTQFMYDEFVEDYEPTKADSYRKKVLDGEEVQIDILDTAGQEDYAAIRD NYFRSGEGFLLVFSITEHESFTATAEFREQILRVKAEEDKIPLLVGNKSDLEERRQVPVEEARSKAEWGVQYVETSAKTRA NVDKVFFDLMRIERTKKMSENKDKNGKKSSKNKKSFKERCCLL

靶点信息

研究背景	This gene encodes a GTP-binding protein that belongs to the small GTPase superfamily and Ras family of proteins. GTP-binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors.
基因ID	5899

基因名	RALB
Swiss	P11234
别名	RALB

产品验证



Western blot analysis of RALB expressed in Mouse lung,Rat lung,SKOV3 using RALB Rabbit pAb at 1:100 0. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn