

PKR/EIF2AK2 Rabbit pAb

货号: B13598

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IF: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:1000 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	57kDa/62kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human PKR/EIF2AK2 (NP_002750.1).
序列	MAGDLSAGFFMEELNTYRQKQGVLKYQELPNSGPPHRRFTFQVIIDGREFPEGEGRSKKEAKNAAAKLAVEILNKEKKA VSPLLLTTTNSSEGLSMGNYIGLINRIAQKKRLTVNYEQCASGVHGPEGFHYKCKMGQKEYSIGTGSTKQEAKQLAAKLAY LQILSEETSVKSDYLSSGSFATTCESQNSLVSTLASESSSEGDFSADTSEINSNSDSLNSSLLMNGLRNNQRKAKRSL APRFDLPDMKETKYTVDKRFGMDFKEI

靶点信息

研究背景	The protein encoded by this gene is a serine/threonine protein kinase that is activated by autophosphorylation after binding to dsRNA. The activated form of the encoded protein can phosphorylate translation initiation factor EIF2S1, which in turn inhibits protein synthesis. This protein is also activated by manganese ions and heparin. Three transcript variants encoding two different isoforms have been found for this gene.
基因ID	5610
基因名	EIF2AK2
Swiss	P19525
别名	EIF2AK1;PKR;PPP1R83;PRKR;Eif2ak2;EIF2AK2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn