

CLTA Rabbit pAb

货号: B13610

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Cynoglossus semilaevis
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	17kDa/23kDa/24kDa/25kDa/27kDa
实测分子量	33kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,NCI-H460,HeLa,A-549
细胞定位	Cytoplasmic side,Cytoplasmic vesicle membrane,Membrane,Peripheral membrane protein,coated pit
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human CLTA (NP_001824.1).
序列	MAELDPFGAPAGAPGGPALGNGVAGAGEEDPAAAFDAQESEIAGIENDEAFAILDGGAPGPQPHGEPPGGPDVAVDGV MNGEYYQESNGPTDSYAAISQVDRLQSEPESIRKWREEQMERLEALDANSRKQEAWEKKAIKELEEYARQDEQLQKT KANNRAAEEAFVNDIDESSPGTEWERVERLCDFNPKSSKQAKDVSRMRSVLISLKQAPLVH

靶点信息

研究背景	Clathrin is a large, soluble protein composed of heavy and light chains. It functions as the main structural component of the lattice-type cytoplasmic face of coated pits and vesicles which entrap specific macromolecules during receptor-mediated endocytosis. This gene encodes one of two clathrin light chain proteins which are believed to function as regulatory elements. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 8 and 12.
------	--

基因ID	1211
基因名	CLTA
Swiss	P09496
别名	CLTA;LCA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn