

TMED9 Rabbit pAb

货号: B13620

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	27kDa
实测分子量	27kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,HeLa,Raji,Mouse liver,Mouse brain,Mouse lung,Rat liver
细胞定位	Endoplasmic reticulum membrane,Endoplasmic reticulum-Golgi intermediate compartment membrane,Golgi apparatus,Single-pass type I membrane protein,cis-Golgi network membrane,trans-Golgi network membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human TMED9 (NP_059980.2).
序列	MAVELGVLLVRPRPGTGLGRVMRTLLLVWLATRGSAFYFHIGETEKKCFIEEIPDETMVIGNYRTQLYDKQREEYQPATPGLGMFVEVKDPEDKVILARQYGSEGRFTFTSHTPGEHQICLHSNSTKFSLFAGGMLRVHLDIQVGEHANDYAEIAAKDKLSELQLRVRQLVEQVEQIQKEQNYQRWREERFRQTSESTNQRVLWWSILQTLILVAIGVWQMRH

靶点信息

研究背景	
基因ID	54732

基因名	TMED9
Swiss	Q9BVK6
别名	TMED9;GMP25;HSGP25L2G;p24a2;p24alpha2;p25

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn