

# AIF1/IBA1 Rabbit pAb

货号: B13639

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>IHC:</b> Mus musculus                                                                             |
| 应用    | <a href="#">IHC</a> <a href="#">IF/ICC</a>                                                           |
| 推荐浓度  | <b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200                                              |
| 理论分子量 | 10kDa/14kDa/16kDa                                                                                    |
| 实测分子量 | 10kDa/14kDa/16kDa                                                                                    |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | THP-1 cells                                                                                          |
| 细胞定位  | Cell projection,Cytoplasm,Cytoplasmic side,Peripheral membrane protein,cytoskeleton,ruffle membrane  |
| 纯化    | Affinity purification                                                                                |

抗原信息

|      |                                                                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human AIF1/IBA1 (P55008).                                    |
| 序列   | MSQTRDLQGGKAFGLLKAQQEERLDEINKQFLDDPKYSSDEDLPSKLEGFKEKYMEFDLNGNGDIDIMSLKRMLEKLGVPKTHLELKKLIGEVSSGSGETFSYPDFLRMMLGKRSAILKMILMYEEKAREKEKPTGPPAKKAISELP |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a protein that binds actin and calcium. This gene is induced by cytokines and interferon and may promote macrophage activation and growth of vascular smooth muscle cells and T-lymphocytes. Polymorphisms in this gene may be associated with systemic sclerosis. Alternative splicing results in multiple transcript variants, but the full-length and coding nature of some of these variants is not certain. |
| 基因ID | 199                                                                                                                                                                                                                                                                                                                                                                                                                                |

|       |                                 |
|-------|---------------------------------|
| 基因名   | AIF1                            |
| Swiss | P55008                          |
| 别名    | AIF1;AIF-1;IBA1;IRT-1;IRT1;lba1 |

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)