

HCN1 Rabbit pAb

货号: B13655

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Rattus norvegicus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	98kDa
实测分子量	99kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,MCF7,OVCAR3,SH-SY5Y,Mouse brain,Mouse eye,Mouse ovary,Mouse lung
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HCN1 (O60741).
序列	MEGGGKPNSSNSRDDGNSVFPKASATGAGPAAAEKRLGTPPGGGGAGAKEHGNSVCFKVDGGGGGGGGGGGGE EPAGGFEDAEGPRRQYGFMQRQFTS

靶点信息

研究背景	The membrane protein encoded by this gene is a hyperpolarization-activated cation channel that contributes to the native pacemaker currents in heart and neurons. The encoded protein can homodimerize or heterodimerize with other pore-forming subunits to form a potassium channel. This channel may act as a receptor for sour tastes.
基因ID	348980

基因名	HCN1
Swiss	O60741
别名	HCN1;BCNG-1;BCNG1;EIEE24;HAC-2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn