

E2F6 Rabbit pAb

货号: B13671

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IF/ICC IP ChIP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 ChIP: 1:50 - 1:200
理论分子量	22kDa/28kDa/31kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-281 of human E2F6 (N P_937987.2).
序列	MSQQRPARKLPSTLLDPTEETVRRRCRDPINVEGLLPKIRINLEDNVQYVSMRKALKVKRPRFDVSLVYLTRKFMDLVRSA PGGILDLNKVATKLGVRKRRVYDITNVLDGIDLVEKKSKNHIRWIGSDLSNFGAVPQQKKLQEELSDLSAMEDALDELIKD CAQQLFELTDDKENERLAYVTYQDIHSIQAFHEQIVIAVKAPAETRLDVPAPREDSITVHIRSTNGPIDVYLCEVEQGQTSNK RSEGVGTSSSESTHPEGPEEEENPQQSEELLEVS

靶点信息

研究背景	This gene encodes a member of a family of transcription factors that play a crucial role in the control of the cell cycle. The protein encoded by this gene lacks the transactivation and tumor suppressor protein association domains found in other family members, and contains a modular suppression domain that functions in the inhibition of transcription. It interacts in a complex with chromatin modifying factors. There are pseudogenes for this gene on chromosomes 22 and X. Alternative splicing results in multiple transcript variants.
基因ID	1876
基因名	E2F6
Swiss	O75461
别名	E2F6;E2F-6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn