

TCOF1 Rabbit pAb

货号: B13688

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	96kDa/144kDa/148kDa/152kDa/155kDa
实测分子量	152kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,293T
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 659-958 of human TCOF1 (NP_001008657.1).
序列	VRVGTQAPRKAGTATSPAGSSPAVAGGTQRP AEDSSSSEESDSEEEKTGLAVTVGQAKSVGKGLQVKAASVPVKGSLGQGTAPVLPGKTGPTVTQVKA EKQEDSESSEESDSEEEAASPAQVKT SVKKTQAKANPAAARAPSAKGTISAPGKVVTAA AQAKQRSPSKVKPPVRNPQNSTVLARGPASVPSVGKAVATAAQAQTGPEEDSGSSEESDSEEEAETLAQVKPSGKTHQ IRAALAPAKESPRKGAAPTTPGKTGPSAAQAGKQDDSGSSSEESDSDGEAPAAV TSAQDQESS

靶点信息

研究背景	This gene encodes a nucleolar protein with a LIS1 homology domain. The protein is involved in ribosomal DNA gene transcription through its interaction with upstream binding factor (UBF). Mutations in this gene have been associated with Treacher Collins syndrome, a disorder which includes abnormal craniofacial development. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	6949
基因名	TCOF1
Swiss	Q13428
别名	TCOF1;MFD1;TCS;TCS1;treacle

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn