

Uhrf2 Rabbit pAb

货号: B13695

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	56kDa/89kDa
实测分子量	95KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 636-733 of mouse Uhrf2 (NP_659122.2).
序列	LQYPAGYPSEKEGKKTGQSKQGSEATKRPASDDECPGDSKVLKASDSTDAVEAFQLTPQQQRLIREDCQNQKLWDE VLASLVEGPNFLKKLEQSFM

靶点信息

研究背景	This gene encodes a nuclear protein which is involved in cell-cycle regulation. The encoded protein is a ubiquitin-ligase capable of ubiquinating PCNP (PEST-containing nuclear protein), and together they may play a role in tumorigenesis. The encoded protein contains an NIRF_N domain, a PHD finger, a set- and ring-associated (SRA) domain, and a RING finger domain and several of these domains have been shown to be essential for the regulation of cell proliferation. This protein may also have a role in intranuclear degradation of polyglutamine aggregates. Alternative splicing results in multiple transcript variants some of which are non-protein coding.
------	---

基因ID	109113
基因名	
Swiss	Q7TMI3
别名	UHRF2;NIRF;RNF107;TDRD23;URF2;Uhrf2

产品验证



Western blot analysis of Uhrf2 expressed in Mouse thymus using Uhrf2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn